


Midi Operating Hysteroscope Set Ø 2.9mm



192mm working length

351-400-043

Operating sheath, with working channel for use with 5 Charr. instruments, with inflow stopcock, for use with outer sheath 351-400-050 and Midi Hysteroscope 351-829-030 or 351-829-130

351-400-050

Outer sheath, Ø 5.0mm, pear-shaped sheath, 192mm working length, with outflow stopcock, for use with operating sheath 351-400-043



351-905-105

Disposable Silicone Valve with LuerLock connection for use with 3 - 10 Charr. instruments sterile packed, package with 10 pieces

EXTENDED working length



235mm working length

351-400-143

Operating sheath, with working channel for use with 5 Charr. instruments, with inflow stopcock, for use with outer sheath 351-400-150 and Midi Hysteroscope 351-829-030 or 351-829-130

351-400-150

Outer sheath, Ø 5.0mm, pear-shaped sheath, 235mm working length, with outflow stopcock, for use with operating sheath 351-400-143